



UCT Specialties, LLC

Technical Data Sheet

Product Name: S1590-KG

Chemical: 3-(N-STYRYLMETHYL-2-AMINOETHYLAMINO)PROPYLTRIMETHOXY SILANE
HYDROCHLORIDE, 40% in MeOH

Cas Number: 34937-00-3

Test	Specifications	Method
Appearance Description	Yellow to golden brown color	SOP82-4.080
Chloride	3.4000 Min.	SOP82-4.069
IR	Pass	SOP82-4.075
Viscosity in cS	1.0000 - 3.0000	SOP82-4.299
Refractive Index N(25,D)	1.3900 - 1.4100	SOP82-4.070
Specific Gravity @ 25 Deg. C	0.8900 - 0.9200	SOP82-4.207

The specification values are only valid in conjunction with the methods listed.

The shelf life for this product is 6 months from the date of shipment.

The shelf life is only valid if the material is stored in the original unopened container.

The specifications and lot analysis are only valid using the listed methods.

The test methods and specifications were current when the product was lotted.