




SpinFiltr[®]

A Whole **New Spin** on dSPE





UCT's new SpinFiltr[®] is taking the hassle out of dSPE

- Simply discard dSPE chamber containing unwanted matrix and sorbent following centrifugation.
- Recover additional sample volume without worrying about the disruption of centrifugation layers.
- Obtain enhanced extract purification with our built-in 0.2 µm filtration device.

Discover a whole new approach to dSPE!

Product Number	Contents
ECQUSF14CT	150 mg MgSO ₄ 50 mg PSA
ECQUSF24CT	150 mg MgSO ₄ 50 mg PSA 50 mg C18
ECQUSF34CT	150 mg MgSO ₄ 25 mg PSA 2.5 mg GCB
ECQUSF44CT	150 mg MgSO ₄ 50 mg PSA 50 mg C18 7.5 mg GCB
ECQUSF54CT	150 mg MgSO ₄ 50 mg PSA 50 mg C18 50 mg Chlorofiltr [®]
ECQUSF64CT	150 mg MgSO ₄ 50 mg C18 50 mg Chlorofiltr [®]
ECQUSF74CT	150 mg MgSO ₄ 50 mg C18

