




A WIDE SELECTION OF INLET LINERS FOR GAS CHROMATOGRAPHS




GC LINERS

GC LINERS


| Confirming Peaks | |
|------------------|---|
| 1 | Pentachloronitrobenzene (IS) |
| 2 | Endrin |
| 3 | 4, 4' DDT |
| 4 | Endrin Aldehyde |
| 5 | Endrin Ketone * 4, 4' DDD * 4, 4' DDE |

* Showed no breakdown from 4, 4' DDT in pesticide breakdown standard.










Before Silanization (Not Deactivated)



After Silanization / Post 50 Injections (Deactivated)



G C LINERS

| Liner Visual | Description / Instrument | Part Number | Quantity |
|--|---|---|-----------------------------|
|  | 2 mm Straight Splitless 2.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCL2MM GCL2MM-5 GCL2MM-25 | 1 Pack 5 Pack 25 Pack |
|  | 2 mm Straight Splitless with Deactivated Glass Wool 2.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCL2MMGW GCL2MMGW-5 GCL2MMGW-25 | 1 Pack 5 Pack 25 Pack |
|  | 2 mm Splitless Gooseneck 2.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCLGN2MM GCLGN2MM-5 GCLGN2MM-25 | 1 Pack 5 Pack 25 Pack |
|  | 2 mm Splitless Gooseneck with Deactivated Glass Wool 2.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCLGN2MMGW GCLGN2MMGW-5 GCLGN2MMGW-25 | 1 Pack 5 Pack 25 Pack |
|  | 4 mm Straight Splitless 4.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCL4MM GCL4MM-5 GCL4MM-25 | 1 Pack 5 Pack 25 Pack |
|  | 4 mm Straight Splitless with Deactivated Glass Wool 4.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCL4MMGW GCL4MMGW-5 GCL4MMGW-25 | 1 Pack 5 Pack 25 Pack |
|  | 4 mm Splitless Recessed Gooseneck with Deactivated Glass Wool 4.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCLRG4MMGW GCLRG4MMGW-5 GCLRG4MMGW-25 | 1 Pack 5 Pack 25 Pack |
|  | 4 mm Splitless Gooseneck 4.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCLGN4MM GCLGN4MM-5 GCLGN4MM-25 | 1 Pack 5 Pack 25 Pack |
|  | 4 mm Splitless Gooseneck with Deactivated Glass Wool 4.0 mm ID x 6.5 mm OD x 78.5 mm Length Instrument: Agilent | GCLGN4MMGW GCLGN4MMGW-5 GCLGN4MMGW-25 | 1 Pack 5 Pack 25 Pack |

DRAWING NOT TO SCALE

PRICES AND TERMS

Our prices are subject to change without notice. The price in effect when we receive your order will apply. All prices are in US Dollars and are F.O.B. Lewistown, PA 17044. Terms of payment are net 30 days.

MINIMUM ORDERS

We welcome all orders, therefore, we do not have a minimum order requirement. When ordering, please include your purchase order number, complete "Ship To" and "Bill To" address, catalog number, quantity, and description of product(s). Also include your name and a phone number where you can be reached should we have any questions concerning your order.

SHIPMENTS

Normal processing is within 24 hours after receipt of an order. Unless special shipping requests have been made, our trained staff will send all orders by UPS Ground service. The appropriate shipping charges (freight & insurance costs) will be added to the invoice, unless otherwise instructed by the customer.

SPECIAL PRICING

We offer special pricing for volume purchases and standing orders. These discounts apply to bonded phase extraction column purchases only. Please call a sales representative for more information on special pricing qualifications.

RETURN POLICY

Our Quality Manager will handle all returns. Before returning merchandise, please call to obtain a return authorization number from the quality manager. We will need to know the reason for the return, date of purchase, purchase order number and invoice number in order to issue a return authorization number. Return merchandise must be received before a credit can be issued. Returns will not be accepted after 90 days. A restocking fee of 25% of the price paid, or a minimum of \$25.00 (whichever is greater) will be charged on all returns.

WARRANTY

All products manufactured by UCT are guaranteed against defects in materials and workmanship for a period of 90 days after shipment. UCT will replace any items that prove to be defective during this time period. The exclusive remedy requires the end user to first advise UCT of the defective product by phone or in writing and must include order number, the lot number and the shipping date.

To initiate this action, photographs of the product, including packaging and labeling of the containers, must be submitted to the UCT Representative for approval. With approval a return authorization can be initiated, and must be received within 30 days. Once the materials arrive at UCT a further inspection of the materials must be completed and accepted by our Quality Manager prior to further action of credits or replacement. UCT's total liability is limited to the replacement cost of UCT products.

This warranty does not apply to damage resulting from misuse.

Placing An Order

Email: info@unitedchem.com

Web: www.unitedchem.com



UCT, Inc.
2731 Bartram Rd.
Bristol, PA 19007

Phone: 800.385.3153
Fax: 215.785.1226
www.unitedchem.com



FORENSICS