




# UCT Enzymes at a Glance



Enzyme Type	Abalonase™	Abalonase™ +	Abalonase™ Ultra	Selectrazyme® (Liquid)	Selectrazyme® (Powder)
<b>Units of Activity</b>	50,000 units/mL	50,000 units/mL	150,000 units/mL	100,000 units/mL	1,000,000-3,200,000 units/gram
<b>Purified</b>	YES	YES	YES	NO	NO
<b>Hydrolysis Buffer Included</b>	YES	YES	YES	NO	NO
<b>Optimal Hydrolysis pH</b>	4.5	5	4.5	5	5
<b>Stability *</b>	1 YEAR	1 YEAR	1 YEAR	1 YEAR	3 YEARS
<b>Storage Conditions</b>	2°C to 8°C	2°C to 8°C	2°C to 8°C	2°C to 8°C	- 20°C
<b>Average Hydrolysis Time **</b>	15-90 minutes	15-90 minutes	0-60 minutes	1-3 hours	1-3 hours
<b>Added Sulfatase Activity</b>	NO	YES	NO	YES-Minimal (Inherent secondary activity)	YES-minimal (Inherent secondary activity)

\* When stored properly-unopened

\*\* Compound dependant



6105-02-01



800.541.0559 [www.unitedchem.com](http://www.unitedchem.com)