



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** D5480-KG

**Chemical:** 1-(DIMETHYLCHLOROSILYL)-2-(P,M CHLOROMETHYLPHENYL)ETHANE

**Cas Number:**

Test	Specifications	Method
Appearance	Pass	SOP82-4.080
IR	Pass	SOP82-4.075
Purity (SUM OF ISOMERS) by GC	94.0000 Min.	SOP82-4.108

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life for this product is months from the date of shipment.*

*The shelf life is only valid if the material is stored in the original unopened container.*

*The specifications and lot analysis are only valid using the listed methods.*

*The test methods and specifications were current when the product was lotted.*